



REO-1 is designed as a light weight, portable independent stand alone instrument. The Graduated Syringe Applicator, are manufactured in accordance with ISO-594 Standards. The Syringe Applicator deposits the fluid sample accurately and consistently on the Ferrogram Substrate or S.E.M. surface at a controlled rate.



ROLLING CONTACT WEAR (tempering problem):
Particles applied to substrate by Syringe Applicator & REO-1 Instrument
(as seen through FERRO-VIEW BX40)

Item:	Description:	Price U.S.D.
Syringe Applicator INST-110-02	High density polymer; Silicon coated to reduce sample cling; For use with the REO-1 Instrument (Sold in 50, 100 or 200 units)	Call for Pricing